

# 恒温水浴锅

## THERMOSTATIC WATER BATH

### 用途特点

恒温水浴锅广泛用于干燥、浓缩、蒸馏、浸渍化学试剂及生物制品，可生物遗传、病毒、水产、环保、医药、卫生、生化试验室、科研教学的必备工具。

- 外壳优质304米粒不锈钢模具冲压而成
- 内胆为优质304不锈钢模具冲压而成
- 保温为整体发泡而成，保温性能极佳

### Characteristics


Thermostatic Water Bath is widely used in drying, concentrating, distilling, dipping chemical reagents and biological products, it is an essential tool for biological genetics, virus, aquatic products, environmental protection, medicine, hygiene, biochemistry laboratory, scientific research and teaching.

- The shell is made of high quality 304 stainless steel mould through stamping
- The inner is made of high quality 304 stainless steel mould through stamping
- The insulation layer is made of the whole foaming and insulation performance is excellent



整体发泡，一次成型  
Overall foaming, one-time molding

### 尺寸图 Dimensions



双孔(2 holes)ZX-S22

### 基本参数 Parameters

参数Parameter 型号Model	温度范围 Temperature range	控制精度 Temperature accuracy	内胆尺寸 Internal size(mm)	外尺寸 External size(mm)	加热功率 Heating power(w)
双孔(2 holes)ZX-S22	室温-99℃	0.5-2℃	300×160×130	345×200×210	800

### 细节图 Details



304不锈钢模具成型  
304 stainless steel molding



304不锈钢隔板  
304 stainless steel plate



304不锈钢加热管  
304 stainless steel pipe